

Key to Action Buttons

Row Actions menu (on rows)

Global Actions menu (top of the screen)

Component Actions menu (on headings)

Additional Quick Reference Guides, training documents, and the Construction and Materials User Guides are available from ODOT.

Role(s): Central Lab Supervisor, Central Lab Manager, Central Lab Admin, and Field Lab Manager.

Accessing the Maintain Test Queue Component:

- 1. From the **Dashboard**, in the **Materials Worksheets** component, click the **Maintain Test Queue** link.
- 2. From this queue, users can view sample record tests from sample records with the following Sample Statuses:

a. Logged

• A Sample Record with this Sample Status will have Tests with a Test Status of **10 - Test in Queue**. Tests have been assigned, but the sample record has not been received.

b. Received at Destination Lab

• A Sample Record with this Sample Status will have Tests in the Receive at Lab Unit queue with a Test Status of **10 - Test in Queue**. The sample record has been received at the Destination Lab but has not been received into a Lab Unit.

c. Received at Lab Unit

• A Sample Record with this Sample Status will have Tests in the Enter Test Results queue with a Test Status of 10 - Test in Queue. The sample record has been received into a Lab Unit.



d. In Testing

A Sample Record with this Sample Status may have Tests in the Enter Test Results queue with Test Statuses of 0 – Test Cancelled, 10 - Test in Queue, 11 – Test Requeued, 20 – Being Retested, 21 – Retest, 40 – Test Complete, or 81 – Test Referred Back to Lab.

e. Pending Authorization

A Sample Record with this Sample Status will have Tests in the Review Samples queue with a Test Status of
60 – Test Approved by Level 2.

Viewing a Sample Record Summary:

1. On the **Maintain Test Queue** component, in the **Sample ID** column, right-click on a link for the option to open the Sample Record Summary in another browser tab.

Viewing a Sample Record Test Summary:

1. On the **Maintain Test Queue** component, in the **Test Number** column, right-click on a link for the option to open the Sample Record Test Summary in another browser tab.



Changing the Indicator if the Test should Count Toward Material Acceptance Actions:

- 1. On the **Maintain Test Queue** component, in the **MAA** column for the desired sample record test, double-click to display the check box.
- 2. Select the check box to indicate **Yes** in the column to consider the test for satisfying material acceptance actions.
 - or -

Clear the check box to indicate **No** in the column to not consider the test for satisfying material acceptance actions.

Changing the Assigned Lab Unit:

- 1. On the Maintain Test Queue component, in the Lab column, double-click to display the check box.
 - **Note:** This will allow you to select an unqualified lab unit for the test.

Changing the Priority:

- 1. On the **Maintain Test Queue** component, in the **Priority** column, double-click to display the choices: High, Medium, and Low.
 - **Note:** This is for information only. There is no system logic that impacts the workflow.



Changing the Planned Test Start Date:

- 1. On the **Maintain Test Queue** component, in the **Planned Test Start Date** column, double-click to type a planned test start date in MM/DD/YYYY format.
 - Note: This is for information only to be added with user discretion. There is no system logic that impacts the workflow.

Changing the Due Date:

- 1. On the **Maintain Test Queue** component, in the **Due Date** column, double-click to type a planned test due date in MM/DD/YYYY format.
 - Note: This is for information only to be added with user discretion. There is no system logic that impacts the workflow.

Adding Test Runs:

- 1. On the Maintain Test Queue component, click the row Actions menu.
- 2. In the **Tasks** section, click the **Add Test Run** task. The system creates additional test runs in the **Enter Test Results** queue, sets the **Test Status** to **10 Test in Queue**, and the **Test Number** field for the additional test runs are numbered as an incremental whole number (for example, 2.0 or 3.0)



Assigning Tests to a Sample Record:

- 1. On the Maintain Test Queue component, click the row Action menu.
- 2. In the **Tasks** section, click the **Assign Tests** task.
- 3. In the **Select Default Tests** section, search for and select any default tests to assign to the sample.
- 4. Click the **Next** arrow in the lower right corner of the screen.
- 5. In the **Select Optional Tests** section, search for and select any optional tests to assign to the sample.
- 6. Click the **Next** arrow in the lower right corner of the screen.
- 7. In the **Add Test Runs and Confirm** section, in the **Test Runs** field, enter the number of test runs to perform for each assigned test on the sample.
- 8. Click the **Confirm** button.
- 9. View the information in the assigned sample record test row.

Splitting a Sample Record:

- 1. On the **Maintain Test Queue** component, on a sample record test row for the desired sample record to split, click the row **Actions** menu.
- 2. In the **Tasks** section, click the **Split Sample Record** task.
- 3. Enter the number of additional sample records you want to create from the split.

Prepared by Info Tech Operating, LLC, DBA Infotech. The AASHTOWare Project logo, the AASHTOWare Project Construction & Materials logo, and the Mobile Tester logo are registered trademarks of AASHTO. AASHTOWare Project, AASHTOWare Project, AASHTOWare Project Construction & Materials, and AASHTOWare Project Mobile Tester are proprietary software products of AASHTO. Mobile Inspector is a proprietary software product of Info Tech Operating, LLC.





- 4. Click Split Sample Record.
- 5. Click **Save**.
 - Note: You can split an unauthorized sample record to create one or more editable copies.
 - Note: When you split a sample record, the system copies information entered on the General, Sample Location, Additional Information, and Sources/Facilities tabs except for the value in the Represented Quantity field. In addition, the system copies sample record remarks, attachments, links, and any contract project item associations. Tests are copied, but test results are not.